

## Method 0040 Sampling Kit #SK-0040 Packing List

<i>SK-0040 includes:</i>				
Part #	Description		Quantity	Qty Shipped
VSB2	VersaCase 2 with panel		1	
PC-3/4E	Extended VersaCase Probe Clamp, Holds MFH-25		1	
MPS-36H	36 inch Heated Probe Assembly		1	
MPL-36G5	36 inch Pyrex Liner, 1/4 inch OD		3	
GF-25T	Filter, PTFE/Glass Fiber		1	
VC40-SET	Clamp Set, M-0040 for VersaCase		1	
<b>VG40-CGS</b>	<b>Glassware Set for Method 0040</b>		<b>1</b>	
<i>Includes:</i>	<b>VG-1/18</b>	3-way Valve with #15 Threads Inlet & purge, 8mm outlet	3	
	<b>VG-34/18</b>	SuperVOST Coil Condenser, 2 inch dia., #15 Joints	1	
	<b>VG-125/18WJ</b>	Knock-out Flask, 125ml with #18 Threads	1	
	<b>VG-2/18</b>	Drying Tube, Charcoal Trap with Glass Grit, #15 Threads	2	
	<b>GA-GL-18B</b>	Bored Cap, #18	13	
	<b>GA-GL-18C</b>	Solid Cap, #18 Threads with Seal	1	
	<b>GA-GL-18T8</b>	TFE Seal Ring, #18/8	3	
	<b>GA-GL-18T6</b>	TFE Seal Ring, #18/6	10	
	<b>TPFA-4/2-047W</b>	TFE Tubing, 1/4 inch OD 1/8 inch ID, per ft.	18	
	<b>TL-7/5</b>	Surgical Tubing, 7/16 inch OD, 5/16 inch ID, per ft.	7	
	<b>4ET4-PFA</b>	Tee, Union, 1/4TU, PFA	1	
MVH-1	3-Way Valve Heater Assembly 50 Watts		1	
SB-3	Insulated Coolant Reservoir		1	
MM5-P	Submersible Coolant Pump, 120V		1	
VU-30	30 Foot VersaCase Umbilical Cable		1	
TBE-6G	Tedlar Bag Enclosure, Uses 25L Bags, 2 Bored SS Unions		1	
TB-18X24-P	25L Tedlar Bag, Poly 2N1 Valve with Septum		3	
MFH-25CA	Heated Filter Assembly, 120V		1	
<i>Includes:</i>	<b>MFH-25C</b>	Clamp, Block, Alum, Threaded	1	
	<b>TFA-25B</b>	Filter, Body, TFE, 25mm, 1/4TU	1	
	<b>MH-50</b>	Heater, Subassy, 50W, w/TC	1	

